

# MICOTOX LTDA



## DETERMINATION OF OCHRATOXIN A (OTA) IN CEREAL GRAINS BY THIN LAYER CHROMATOGRAPHY (TLC)

AOAC Official Method 973.37

### EXTRACTION

1. Weigh out 25 g of ground sample into a 250 mL Erlenmeyer flask.
2. Add 12.5 mL of phosphoric acid (H<sub>3</sub>PO<sub>4</sub>) 0.1M (3.42 mL of 85% phosphoric acid in 500 mL water) and swirl to wet sample.
3. Add 125 of methylene chloride (CH<sub>2</sub>Cl<sub>2</sub>) and shake for one hour on gyratory shaker.
4. Filter through quantitative filter paper (e.g. S&S 589-3) and pipet exactly 5 mL of extract into a 15 x 85 mm test tube.

### PURIFICATION

5. Place a solid phase Micotox® M2200 ochratoxin A column onto a solid phase extraction apparatus. Pour the 5 mL extract into the column (no preconditioning is necessary). Drain by gravity (<2.0 mL/min)\* and discard solution.

\*NOTE: Solution has to drain through the column at a maximum flow rate of 2 mL/min. Before using the M2200 Micotox® column, press firmly the top frit of the column with a plunger. This will prevent the sample extract from draining through the column too fast.

6. When top of solution reaches top of packing add 5 mL of hexane followed by 15 mL of CH<sub>2</sub>Cl<sub>2</sub>, added in three successive volumes of 5 mL each. Discard solution.
7. Elute the ochratoxin A with 20 mL of CH<sub>2</sub>Cl<sub>2</sub>/HCOOH (formic acid) 99+1 in two successive volumes of 10 mL. Use an appropriate test tube or container to collect the eluate (e.g. a 40 mL capacity conical test tube).
8. Evaporate the eluate to dryness under nitrogen or in a 60°C waterbath using vacuum.

### THIN LAYER CHROMATOGRAPHY

- Dissolve the dry residue with 100 µL toluene-acetic acid, 99+1 (add the solvent mixture, stopper the test tube and mix well on vortex for 30 seconds).
- Spot 20 µL of each sample along with 5, 10, and 20 µL of working standard on a 10 x 10 cm silicagel 60 plate (Merck 1.05553). Working standard contains 1.0 µg/mL ochratoxin A in toluene-acetic acid, 99+1 (available from Micotox Ltda.).
- Develop plate in toluene-methanol-acetic acid (90+5+5) to about 1 cm from the top of the plate. Let plate air dry in a fume hood.
- View TLC plate under long wave UV light (365 nm) and determine ochratoxin A content by comparison with the standard spots. Ochratoxin A has a blue-green fluorescence and its approximate R<sub>f</sub> is 0.6.

### CALCULATIONS

25 g sample are extracted with 125 mL extraction solvent; 5 mL from the extract are loaded onto column, taken to dryness and then dissolved with 100 µL toluene-acetonitrile; 20 µL are spotted on TLC plate.

Sample equivalent in g is:

$$25 \text{ g} \times 5/125 \text{ mL} \times 0.02/0.1 \text{ mL} = 0.2 \text{ g}$$

$$\text{ng/g (ppb)} = \frac{\text{ng of toxin on plate}}{0.2 \text{ g}}$$

Amount of ochratoxin A standard spotted on plate (ng):

Vol. (µL)	Ochratoxin A
5	5 ng
10	10 ng
20	20 ng

Equivalent concentration of ochratoxin A standard on TLC plate (ng/g):

Vol. (µL)	Ochratoxin A
5	25 ppb
10	50 ppb
20	100 ppb

Limit of detection: <5 ppb.

Limit of quantitation: 10 ppb

Confirmation: spray TLC plate with an alcoholic solution of sodium bicarbonate (6.0 g of NaHCO<sub>3</sub> in 100 mL H<sub>2</sub>O + 20 mL ethanol) and let dry. View under long wave UV light. The blue-green fluorescence of ochratoxin A turns to blue and becomes more intense.

OCHRATOXIN A by TLC - MICOTOX